



Excel Engineering Co., LTD.
Yobe - ru Co., LTD.

General Catalog



Contact us:

EE-2019-8

Excel Engineering Co., LTD.

Seika Building 9F, 2-12-14 Ogawamachi, Chiyoda-ku, Tokyo, Japan 101-0052
Phone: 03-5280-7120 Fax: 03-5280-7123
<http://www.excel-jpn.com>

Yobe-ru Co., LTD.

Seika Building 9F, 2-12-14 Ogawamachi, Chiyoda-ku, Tokyo, Japan 101-0052
Phone: 03-5280-7120 Fax: 03-5280-7123
<http://www.yobe-ru.jp>

- For wireless homecare solution products published on this catalog, it is a premise that only one system is being used in a house
- For the products that use the radio wave, the detection signal reaching ranges may vary based on an environment. As a result of that, the ranges may become short or the devices may not function at all. For the products that use a 2.4 GHz frequency band, they are prone to develop a radio signal interference since they use the same signal frequency band with other devices such as a wireless LAN device.

- Colors of the product pictures on this catalog may be different due to a printing process
- Design and specifications are subject to change without notice

Care Communication

Excel Engineering helps provide independent living for care recipients and an effective, efficient and caring environment for caregivers. Our products help support a streamlined communication environment for all involved.

Our hope is that both care recipients and caregivers have meaningful and enjoyable lives together.



Each page has a QR code. When you scan it you will see an instructional manual for the product.

Table of Contents

Total Care Solutions for Senior Care Facilities

| | | | |
|---|-------|-------|--|
| iBaby Senior · iBaby Care | pages | 3-4 | |
| Wireless Nurse Call System: NC-9000 | | 5-6 | Wireless Nurse Call System (920MHz Band/LoRa) |
| Wireless Nurse Call System: NC-5000 | | 7-8 | Wireless Nurse Call System (2.4GHz Band) |
| Cloud-Based Total Care System | | 9-10 | Cloud-Based Total Care System |
| Detection Sensor's List | | 11-12 | |
| Floor Sensor | | 13 | Wireless, cordless, compact sensor for wandering prevention |
| Bed Sensor | | 14 | Help prevent for falling accidents |
| Wandering/Leaving-Bed Plate Detection Sensor: EPS-50 | | 15 | Early detection and alert system for a potential accident |
| Side-Bed Sensor Pad | | 16 | |
| Side-Rail Sensor Pad | | 16 | Install on a side rail of a bed |
| Side Bed Sensor Pad | | 16 | |
| Magnet Sensor | | 16 | |
| Infrared Double Sensor Detection Monitor System: MS-3100D | | 17-18 | Infrared Detection Monitor System for wandering/falling prevention |
| Wheelchair Sensor Pad | | 19 | Pad sensor for wheelchairs |
| Magnet Sensor | | 20 | Magnet sensor for wheelchairs |
| Pad-less Wheelchair Sensor | | 20 | Pad-less sensor |
| Toilet Sensor | | 21 | Help detecting falling accidents in the bathroom |
| Diaper Sensor | | 21 | Help detecting when the sensor gets wet |
| Door Sensor | | 22 | Detect when a door opens |
| Wireless Nurse Call Button | | 23 | |
| Wired Nurse Call Activation/In-room Alert System: EXC-02B | | 23 | |
| Accessory Options | | 24 | Accessory products for the detection sensors |
| Adapter Box | | 25-26 | Adapter box that can be connected with existing nurse call systems |

Homecare Solutions

| | | |
|---|-------|---|
| *Indicates long-term care insurance approved products for rent. | | |
| SRH9000* | 27-28 | Long range wireless system that does not require repeaters |
| SH6900* | 29-30 | Live video monitoring/infrared detection sensor |
| Homecare Wireless Nurse Call System: NRM-6900 | 31 | Two-way voice communication/live video monitoring system |
| Homecare Wireless Nurse Call System: HT-I50 | 32 | Ornament-style transmitter can alert caregiver by sound and light |

Specifications and Other

| | |
|-----------------------|-------|
| Specifications | 33-35 |
| Company's Information | 36 |

The price is FOB Japan for 1 pic

NEW iBaby Senior



Single item & External Sensor Package Prices



iBaby Senior
RH7100
TAIS:00682 - 000078

\$210



iBaby Senior
with Floor Sensor
RH7100
TAIS:00682 - 000079

\$285



iBaby Senior with Bed
Sensor
RH7100
TAIS:00682 - 000080

\$270



iBaby Senior with Plate
Sensor
RH7100
TAIS:00682 - 000081

\$270



iBaby Senior with
Infrared Sensor
RH7100
TAIS:00682 - 000082

\$385

NEW iBaby Care




Single item & External Sensor Package Prices



iBaby Care
RH7000

\$210



iBaby Care
with Floor Sensor
SH7000

\$285



iBaby Care with Bed
Sensor
SR7000

\$270



iBaby Care with Plate
Sensor
SP7000

\$270



iBaby Care with
Infrared Sensor
SK7000


\$385

NC-9000 Wireless Nurse Call System


The wireless signal range could reach more than one kilometer.

Home Environment Usage Examples


1 The system is suited to use in a multi-story high raising building



2 The system is suited to use between buildings



3 The system is suited to use within a neighborhood. When the system is used in a residential area the wireless signal can reach up to several hundred meters.




*Because of the technology used in this system, the wireless signal range varies based on the network environment. Please test the wireless signal strength and range before you use the system.

Facility Usage Examples

- Group Homes
- Small Multifunctional Facilities
- Assisted Living Facilities
- Retirement homes
- Geriatric Health Care Facilities, etc.

Without using repeaters or signal extenders, the transmitter and the receiver can be directly connected wirelessly. This will help alleviate the financial burden for a facility.



Sensors:
Please refer to
page 11 and 12.

Nurse Call Button
NC-9000TXPB

\$ 225

Transmitter
NC-9000TXBF

\$ 225

NC-5000:

- Sensors
- Alert button

Signal Frequency
Converter
NC-9000GW

\$ 500

Portable Receiver
NC-9000PRX

\$ 260

- Can register up to 100 channels
- Rechargeable

Base Station Receiver
NC-9000RX

\$ 270


- Can register up to 16 channels
- AC adapter

Alert Display
Monitor/Receiver
NC-9000DRX


\$ 1590

- Can register up to 521 channels
- Automatically save reports and alerts
- Can manage time periods and users
- Select and filter contents of reports and print them

例



〇月〇日~〇月〇日
Report History



〇月〇日~〇月〇日
Alert Report History


Wireless Nurse Call System NC-5000

Wireless Nurse Call System that adapts various facility needs!

① Excel Engineering' s various detection sensors can be connected to the existing nurse call system wirelessly.

② Detection and/or nurse call request signals can be transmitted to the base and/or the portable receivers.

NC-5000:Installation Example



Wireless Nurse Call System NC-5000 Products

Wireless Adapter Box

(NC-5000NRX-10~22) for an existing nurse call



\$ 150



- This wireless adapter box has capabilities to receive detection signals and connect to an existing nurse call system. The wireless environment will help in creating a safe (no tripping over the cords) and more accessible caring environment.
- It is used in facilities which have a nurse call system
- Up to 3 transmitters can be registered.
- An alarm sound can be stopped temporarily.
- When the transmitter's battery become low the adapter box's LED light starts flashing when it receives a detection signal. Please check this at least once a month).
- Please refer to page 21 for the adapter box's connector

Base Station Receiver

NC-5000RX(C)



\$ 240



- A caregiver can be alerted by "alert" lights, sound, or melodies (8 tunes)
- Up to 16 different detection sensors and/or communication devices can be connected to this receiver
- Multiple receivers can be used simultaneously
- When the transmitter's battery becomes low the LED light starts flashing when an alert signal is received. Please confirm if the system is working properly at least once a month.

Portable Receiver

NC-5000PRX2



\$ 225



- A caregiver can be alerted by an "alert" sound or vibration
- Up to 63 different detection sensors and/or communication devices can be connected to this portable receiver
- Multiple portable receivers can be used simultaneously
- When the transmitter's battery becomes low the LED light starts flashing when an alert signal is received.

Wireless Repeater

NC-5000RP(C)



\$ 150

- This repeater increases the wireless network coverage area
- By connecting multiple repeaters, it is possible to cover the wireless network in an entire facility
- The wireless coverage can be identified by pressing the test button

Mini Transmitter

NC-5000TXBF




\$ 90

This mini transmitter can be used for the bed sensor or the floor sensor. The delayed alert timer can be set at 0, 1, 2, or 4 seconds

Portable "nurse call" button for the wireless

adaptor: NC - 5000TXPB



\$ 110

This is a small portable "nurse call" button

Portable Pause Button Box for the wireless

adaptor: NC - 5000TXPH



\$ 100

This is a portable pause button box for the wireless adaptor. When a caregiver wants to clean the room and/or make the bed for a care recipient a nurse call system can be temporarily deactivated by pressing the pause button. The nurse call system can be reactivated by pressing the button on the device. By using the device, false activation of the nurse call system can be prevented. After 3 minutes elapse, the "pause" function is deactivated automatically. This function can prevent someone from forget to deactivate the "pause" function

Mini Repeater



NC-5000RP2



\$ 225

- This repeater increases the wireless network coverage area
- By connecting multiple repeaters, it is possible to cover the wireless network in an entire facility
- To use the repeater, Simply plug it in to an electrical outlet

Detections sensors that support NC-5000 system



*To see more detailed information of each sensor, please refer to the corresponding pages

Since the NC-5000 is a wireless signal transmission system, we will examine your wireless signal environment before preparing an estimate and installation of the system. Based on your facility's environment conditions, the system may require a repeater(s). Please examine the wireless signal coverage in your facility. Please contact us for more detailed information on the products. We are looking forward to hearing from you.

Cloud-Based Total Care System

A care recipient can be monitored by using various wireless detection sensors on your PC or smartphone within your network.

- Cost effective installation! As long as you have a Wi-Fi environment is available in your home the initial cost is minimal. You can rent our detection sensors you can keep the budget lower.
- The Vital monitoring sensor which has a high demand in a senior care facility can be used with this system.
- Easy data sharing! Care recipients' detection information can be shared between a care facility and the care recipient's family members.

PC and Smartphone



Internet



Sensor Devices

Receiver (Built in Router)

- It is compact and can be installed in a very small space. The wireless signals can be detected within the radius of 15m.
- By installing multiple receivers, a care recipient's position within a facility can be identified.



\$315

Pulse Sensor

- This sensor can measure an individual's pulse count
- The built-in tri-axial accelerometer can detect an individual's motions
- There is a wireless "push button" built in to alert a caregiver

\$315



NEW Twin Sensor Package

- Two detection sensors can be combined to detect "getting up" from the bed, "Leaving" the bed, and presence of an individual in bed.
- The prices vary based on the combination of the sensors

Leaving-Bed Detection Sensor

- There are two sensors available: Bed Sensor and Floor Sensor
- The sensors detect when an individual is trying to leave the bed



\$225

Diaper Sensor

- The sensor alerts when the sensor gets wet.
- This will help a caregiver to know when to change the diaper



\$200

Nurse Call Button (Wrist Band Type)

- The built-in tri-axial accelerometer can detect an individual's motions
- There is a wireless "push button" built in to alert a caregiver
- This device is IPX7 water proof rating product

\$200



Wandering/Leaving-Bed Plate Detection Sensor


- This device is very light and compact design
- The sensor detects wandering / leaving-bed incidents

\$295




[Sold separately: Sensor Tate Built-in Pads (30 counts)]









System Schema



Caring App for your smartphone



NEW Twin Sensor Package Usage Option


| | Usage Options | Sensor activation status and a care recipient's position status | | | | A care recipient's position status | |
|----------|--|---|---|---|---|------------------------------------|--|
| Option 1 |  2 Plate Detection sensors + 1 Floor sensor | Sensor A | ✓ | ✓ | ✗ | ✗ | Position Status 1  |
| | | Sensor B・C | ✗ | ✓ | ✓ | ✗ | |
| | | Position Status | 2 | 2 | 4 | 5 | |
| Option 2 |  2 Plate Detection sensors | Sensor A | ✓ | ✓ | ✗ | ✗ | Position Status 2  Position Status 3  |
| | | Sensor B | ✗ | ✓ | ✓ | ✗ | |
| | | Position Status | 2 | 2 | 3 | 5 | |
| Option 3 |  1 Plate Detection sensor + 1 Floor sensor | Sensor A | ✓ | ✓ | ✗ | ✗ | Position Status 4  Position Status 5  |
| | | Sensor C | ✗ | ✓ | ✓ | ✗ | |
| | | Position Status | 3 | 4 | 1 | 1 | |

- A caregiver can monitor a care recipient via a PC or a smartphone 24/7
- A large number of care recipients' caring data can be managed.
- There are various alert settings available (e.g., upper and lower limits of pulse counts, an emergency call request, duration of inactive status)
- Based on the location of the receiver, a care recipient's position within a care facility or home can be identified.
- An alert signal is sent to your PC or smartphone along with corresponding information from the various detection sensors

- A 3-year customer support fee is included and the cloud service is free
- Within the 3 year of contract, if any of the devices become malfunctioning we will replace them for free
- We do have demo devices. You are more than welcome to try out our product
- A PC and/or a smartphone are not included

List of Sensors


Around the Bed



Floor Sensor

Can work with NC5000 and NC9000

P13



Bed Sensor

Can work with NC5000 and NC9000

P14





Plate Sensor

Can work with NC5000 and NC9000


P15



Infrared Sensor

Can work with NC5000 and NC9000


P17



Side-Bed Sensor

Can work with NC5000


P16



Side-Rail Sensor

Can work with NC5000 and NC9000

P16




Magnet Sensor

Can work with NC5000 and NC9000


P16

Wheelchair




Wheelchair Sensor Pad

P19



Pad-Less Wheelchair Detection Sensor


P20



Magnet Sensor

P20

Toilet Seat Sensor




Toilet Seat Sensor

Can work with NC5000 and NC9000

P21

Diaper Sensor




Diaper Sensor

Can work with NC5000

P21


Exit and Entrance



Door Sensor

Can work with NC5000

P22






Infrared Sensor

Can work with NC5000

P22

Floor Sensor: For wandering prevention!


- When an individual steps on the sensor pad it detects!
- This sensor pad alerts when an individual's body weight creates a pressure on the pad by the bed (where the individual places his or her feet when he or she leaves the bed).
 - The pad would be most appropriate for individuals who need to be monitored for the potential risk of falling or wandering.



| | Floor sensor only | Wired Floor sensor & a Nurse Call Connection Package The package includes the Hybrid Floor Sensor, an Adapter Box, and Nurse Call connection cables | Wireless Floor sensor & NC5000 Package The package includes the Hybrid Floor Sensor, a wireless Adapter Box, and a Mini Transmitter |
|----------------------------------|---------------------|--|--|
| Hybrid Floor Sensor: Size Small | EFM-57/S-C \$110 | EFMS-57/S \$160 | SNC5000F-57/S \$295 |
| Hybrid Floor Sensor: Size Medium | EFM-57/M-C \$175 | EFMS-57/M \$220 | SNC5000F-57/M \$340 |
| Hybrid Floor Sensor: Size Large | EFM-57/L-C \$200 | EFMS-57/L \$245 | SNC5000F-57/L \$180 |

Bed Sensor: For help prevent an individual from falling off the bed.


- This sensor is designed to help prevent an individual from falling off the bed.
- Based on an individual's needs a location of the sensor is selected. When an individual attempts to get up or the body moves too close to the edge of the bed the sensor alerts.



Bed Sensor Installation Example

- When an individual attempts to get up it detects!
- Place the sensor just underneath of the individual's shoulder blades.
 - When pressure from an individual's body weight decreases a caregiver is alerted

Wired Bed Sensor




EXC-02B is used to connect with an existing nurse call system

Size:Wide EPB54-W-I **\$60**

EXC-02B & Bed Sensor (Size: Wide)
Package:The Bed Sensor can be folded in half.
EXC02S-5W **\$275**

The package includes one Bed Sensor (size Wide), one EXC-02B, One Adapter Box, and Nurse Call connection cables

Wireless Bed Sensor



The Wireless Adapter Box is used to connect a nurse call system

EPB58-W-2 **\$100**


NC-5000 & Wireless Bed Sensor
Package:The Bed Sensor can be folded in half.
SNC5000B-58W-2

The package includes one Wireless Bed Sensor (build-in transmitter) and one Wireless Adapter Box **\$295**

Affordable Sensor Pad



When an individual attempts to get up it detects!



| Product # | Types of Plugs | Dimensions |
|-------------|----------------|------------------|
| EPB90-W | Modular Plug | 381 x 762 x 6 mm |
| EPB90-W/MP | Mini Plug | |
| EPB90-W/RCA | RCA Plug | |

*The duration of warranty on these products is 90 days from purchase; however, they can be used longer. They are not intended to be used for multiple care recipients.

Floor Sensor: For wandering prevention!

The sensor uses the Seesaw Switch Technology which helps to minimize false detections. It is very compact yet, it can detect presence of a care recipient on the bed sensitively without complicated detection setting adjustments. The sensor is durable and invisible when installed. Further, the sensor can be used as a bed sensor and a side-bed sensor. The sensor can be connected to Excel Engineering's Cloud-Based Total Care System, (see page 9-10), NC-9000 Wireless Nurse Call System (see page 5-6), NC-5000 (see page 7-8), and EXC02S Wired Nurse Call Activation/In-Room Alert System (see page 23)




Early detection and alert:

It detects when a body (the majority of the body weight) gets too near the edge of the bed for preventing a potential falling accident



Relationship between the sensor activations and what actually happened to a care recipient at a care facility



440 × 50 × 16 mm

Wandering/Leaving-Bed Plate Detection Sensor:
EPS-50

\$ 65

Wired Nurse Call Connection Package EPSS-50

Package includes:
- One Wandering/Leaving-Bed Plate Detection Sensor
- One adaptor box
- One adaptor cable

\$ 130

Wireless Nurse Call System Package SNC5000PS

Package includes:
- One Wandering/Leaving-Bed Plate Detection Sensor
- One wireless adaptor box
- One NC5000 transmitter

\$ 270


Wired Nurse Call Activation/In-Room Alert System Package EXC02S-PS

Package includes:
- One Wandering/Leaving-Bed Plate Detection Sensor
- One EXC02S unit
- One adaptor box
- One adaptor cable

\$ 280

Side-Bed Sensor Pad

Side-Bed Sensor Pad (Wireless)



The Side-Bed Sensor Pad connected to a nurse call system wirelessly.
*A cord from the sensor can be shortened as an option.


EPB54 \$ 40

Side-Bed Sensor Pad and NC5000 Package: SNC5000F-5BS

\$ 250


The Package includes:
One Side-Bed Sensor Pad, one NC5000 transmitter, and one Adapter Box

Side-Rail Sensor Pad




Side-Rail Sensor Pad - This sensor can be placed over the bed-side-rail. When an individual grabs the side rail (i.e., when the individual creates the pressure on the sensor pad) it activates a nurse call system. Since the transmitter is built inside of the pad, it creates a safe cordless environment around a care recipient. The side-rail sensor pad can be used with Excel's wireless nurse call system, NC5000 (see page 7-8).

When an individual grabs the side it activates a nurse call system.




Cordless Side-Rail Sensor Pad (EPB58-SR/L)



EPB58-SR/L \$ 105

Small Cordless Side-Rail Sensor Pad (EPB58-SR/S)





EPB58-SR/S \$ 50

Magnet sensor

Magnet sensor and NC5000 Wireless Nurse Call System package (SN5000TXMG)

An example of when the Magnet Sensor is used for a



The package includes:
One Magnet Sensor and one NC5000

Wireless

SNC5000TXMG

\$ 170

The Infrared Double Sensor Detection Monitor System (MS3100D)



See the instructional manual




The infrared sensors detect movements on the bed. Two different kinds of sensors (Double-triggered sensor system) detect an individual's departure from the bed, absence from the bed, and potentially hazardous locations on the bed.

MS-3100D \$295



Type A: This setting is used for detecting the movement of an individual on the bed who is "getting up from the bed." Using the two types of infrared sensors helps minimize false alarms!

Two types of infrared detection sensors are used effectively to detect an individual on the bed who is "getting up from the bed".




Type B: This sensor detects when the body gets too near the edge of the bed (antifalling detection). The sensor also detects subtle sluggish movements!

When only the Body-position detection sensor is used the infrared beams stop, and the programmable timer starts (ranges from 0-60 minutes). When the timer stops the infrared beams start trying to detect the presence of the individual on the bed. If the individual does not come back to the detection range (on the bed) in the programed time an alert system will be activated.




Other applications of using the mS3100D

Type C: This setting is used for detecting an individual's potential wandering ("not back in bed") and emergency situations. The programmable timer can be set for 3, 6, 10, 20, 30, or 60 minutes.



Type D: The sensor can be placed by a bed to be used for monitoring if an individual leaves his/her bed.




MS-3100B Accessories



Long stand MS-3000-2

Includes: a long pole, and a stand

\$20



Short stand MS-3000-3

Includes: a short pole, and a stand

\$20



Bed Holder BH-2

Easy to install!

\$15



Adjustable Angle Connector PT-1


Easy to change the angle. It can be attached to the holder or the pole.

\$35

MS-3100D & Nurse call Connector Package (Wired)

MS-3100D is connected to a Nurse Call system

Wired Type




| Package Type | Product # | Contents | |
|------------------------------|------------|---|--------------|
| MC-3100 with the Long Stand | MS3100D-LB | ● MC-3100 (1) ● Long Stand (1) ● 2PNC cable (1) ● Adapter Box (1) | \$340 |
| MC-3100 with the Short Stand | MS3100D-SB | ● MC-3100 (1) ● Short Stand (1) ● 2PNC cable (1) ● Adapter Box (1) | \$340 |
| MC-3100 with the Bed Holder | MS3100D-H | ● MC-3100 (1) ● Bed Holder (1) ● 2PNC cable (1) ● Adapter Box (1) | \$340 |

MS-3100D & Nurse call Connector Package (Wireless)

The Wireless Adapter Box is used to connect to a nurse call system


Wireless Type



| Package Type | Product # | Contents | |
|------------------------------|------------------|--|--------------|
| MC-3100 with the Long Stand | SNC5000-MS3100EL | ● MC-3100 (1) ● Long Stand (1) ● Wireless Adapter Box (1) | \$475 |
| MC-3100 with the Short Stand | SNC5000-MS3100ES | ● MC-3100 (1) ● Short Stand (1) ● Wireless Adapter Box (1) | \$475 |
| MC-3100 with the Bed Holder | SNC5000-MS3100EH | ● MC-3100 (1) ● Bed Holder (1) ● Wireless Adapter Box (1) | \$475 |

Wheelchair Sensor Pad

Falling accidents typically occur when an individual tries to get up from a wheelchair. This sensor pad is designed to use with a wheelchair.



When an individual tries to get up from a wheelchair it detects! Falling accidents typically happen when an individual tries to get up from a wheelchair. When the weight of the individual becomes less on the wheelchair due to the attempt of the departure from the wheelchair the device alerts the individual as well as a caregiver.

Wheelchair Sensor Pad (EPC-53) Anti-bacterial Water resistant

EPC-53 \$ 45



305 × 382 × 6.5mm

Affordable Wheelchair Sensor Pad (EPC-90)

EPC-90




305 × 381 × 6mm

*The duration of warranty on these products is 90 days from purchase; however, they can be used longer. They are not intended to be used for multiple care recipients.

Wheelchair Sensor Pad (EPC-53) & NC5000 Package

By using the wireless Adapter Box, the Wheelchair Sensor Pad can be connected to the nurse call system wirelessly.

SNC5000C-5C \$ 250




Package contents

- One Wheelchair Sensor Pad (EPC-53)
- One NC5000 Transmitter
- One Wireless Adapter Box

Wheelchair Sensor Pad (EPC-53) & EXC02-5C Package

The Wheelchair Sensor Pad can be connected to EXC02-5C. The EXC02-5C can alert an individual in a wheelchair by a melody/buzzer sound when he/she tries to get up. The melody/buzzer sound stops when the individual sits back down.


EXC02S-5C \$ 200



Package contents

- One Wheelchair Sensor Pad (EPC-53)
- One EXC02-5C

Pad-Less Wheelchair Detection Sensor (ETL08-I)



Install the ETL08-I underneath of the wheelchair seat. The ETL08-I can alert an individual in a wheelchair by an alarm/chime/buzzer sound when he/she tries to get up. The melody/buzzer sound stops when the individual sits back down. Less noticeable from the individual in the wheelchair.

Voice Setting:

The voice recording can be up to 15 seconds using your own voice or a favorite comforting voice.

"Please be seated!"

"It is dangerous to get up!"



ETL08-I

\$ 190


Magnet Sensor

Magnet Sensor



When a magnet attached to the main unit (a bed or a wheelchair) separates from the main unit an alert signal is sent wirelessly to a nurse call station or an alarm sound is generated to warn an individual in a wheelchair. Please see the diagrams below.


An example of when the Magnet Sensor is used for a wheelchair.



Magnet sensor and NC-5000 Wireless Nurse Call System package (SN5000TXMG)

Excel's wireless adaptor box is used to connect the Magnet Sensor wirelessly. The detection signal will be transmitted to a nurse call station.

SNC5000TXMG \$ 170




This package includes one Magnet Sensor and one Wireless Adaptor Box

Magnet Sensor only (EMG-01)

An alarm sound warns the individual in the wheelchair. When the magnet is replaced on the main unit the alarm sound will stop.



EMG-01 \$ 50



Toilet Sensor & Disposable Diaper Sensor

Toilet Sensor (ETP56-5)

The sensor pads alert when an individual decreases the body weight pressure from the toilet seat or pad.



The sensor pad is placed underneath of the toilet seat.

To function as an on-site alert system for monitoring an individual's departure from the toilet seat, use the Toilet Sensor and EXC-02B together (see page 12).
To Connect an existing Nurse Call System wirelessly, use the Toilet sensor and NC-5000 (see page 7-8) together.




Water-Resistant
& Antibacterial

Toilet Sensor (ETP56-5)

\$55

Disposable Diaper Sensor

The sensor alerts when the sensor gets wet! This will help a caregiver to know when to change the diaper. This sensor can be connected to an existing nurse call system wirelessly via NC-5000 (see page 7-8).



Sensor Tapes (WCT): The sensor tape is placed inside of the diaper. The tape consists of three layers (garment/aluminum/garment). This makes a user more comfortable when he or she wears the diaper.

65 x 410mm



This is a one-time-use disposable tape.


100 pieces WCT100

400 pieces WCT400

\$10/100pcs
\$45/400pcs

Sensor Tape Built-in Pads (WCP): This built-in sensor tape pad absorbs moisture. There are 30 pads in this item.

210 x 490mm
Absorbable capacity : 500ml



This is a one-time-use disposable tape.

WCP500

\$10

Diaper Sensor Transmitter (NC-5000TXWA): The transmitter can send a detection signal wirelessly to a receiver.



NC-5000TXWA


\$135

Door/Corridor Sensor


When the door opens it detect!

The door sensor/transmitter and the magnet unit should be installed to an individual's room door (see a diagram above). When the individual opens the door the door sensor/transmitter units are forced to separate from each other. This triggers an alert to a caregiver. The Excel's NC-5000 Wireless Nurse Call system can be connected with the door sensor and the NC-5000's various receivers can be activated when the door is opened. It is intended to prevent wandering and to monitor a door's open/shut activities.

Door Sensor



Detect!



Door sensor/transmitter's dimensions:
41mm x 69mm x 9mm
Magnet's dimensions:
12mm x 50mm x 7mm


Door sensor/transmitter NC-5000TXDR

\$145

Includes: One Door sensor/transmitter & one Magnet Unit

Exit and Entrance: Infrared Sensor(RH-5000TXIR)

The infrared Sensor can detect an infrared generated from a person and alert a caregiver via NC5000 (see page 7-8).




38 x 109 x 69mm

Door Sensor RH-5000TXIR

\$145


Portable Wireless Nurse Call Activator: NC-5000TXNCB



This is a portable nurse call activator. Pressing any of the 4 buttons will activate a nurse call system wirelessly. For an individual who has difficulty pushing a regular nurse call button, he or she can squeeze or place his or her hand or elbow on the nurse call activator to activate the nurse call system.

NC-5000TXNCB \$110

- This portable nurse call activator can be connected with Wireless Nurse "Call" System NC-5000 (Refer to page 7- 8).
- Pressing any of the 4 buttons will activate a nurse call system wirelessly
- This device includes a lanyard to carry the activator
- This device can be used for someone who has difficulty pushing a regular nurse call button




Wired Nurse Call Activation/In-room Alert System (EXC-02B)




EXC-02B \$170

To function as the in-room alert system with the Bed Sensor



To function as the wheelchair alert system



The Bed Hanger can be used to hang the EXC-02B on the back of a wheelchair.

This unit can be interfaced with the Excel Engineering's Bed Sensor and activate a nurse call system. In addition, this unit, with the Bed Sensor, can be used as an in-room alert system for facilities without a nurse call system. A timer can be set (i.e., 0,1,2, and 4sec) before the bed sensor pad alerts a nurse call system. This can minimize a false alarm caused by an individual departs from the sensor pad and returns onto the sensor pad momentarily.

The Wheelchair Sensor Pad can be connected to EXC02-5C. The EXC02-5C can alert an individual in a wheelchair by a melody/buzzer sound when he/she tries to get up. The melody/buzzer sound stops when the individual sits back down. There are multiple melodies and sounds available. If multiple users use the system simultaneously they can be identified by the different melodies and sound they use.

Optional Accessories

We have a variety of nurse call connector cables and adaptors to meet your needs.

Nurse call connector cable for MS-3100B: MS-03

\$30

This connector cable can be connected to MS-3100B and an Excel's Adaptor box.



Extension cable for 2P metal connector: CC06

\$30

This an extension cable for 2P metal connector. The normal length of the cable is 2m; however, the length can be customized for a longer length.



Modular jack to 3.5 mm mini plug adaptor: CC09

\$15

This adaptor can be connected to a detection sensor with a modular jack and to a 3.5 mm mini plug.



2P metal connector to 3P metal connector adaptor: CC11

\$30

This adaptor can be connected to a detection sensor with a 2P metal connector and to a 3P metal connector.



Modular jack to 2P metal connector conversion adaptor: CC14

\$20

This adaptor can be connected to a detection sensor with a modular jack and to a 2P metal connector.



3.5 mini plug to 2P metal connector conversion adaptor: CC16

\$25

This adaptor can be connected to a detection sensor with 3.5 mini plug and to a 2P metal connector.



Hybrid Sensor cable for Modular jack CC57-1

\$15

This cable can be connected to a hybrid detection sensor and to a modular jack.



Nurse call connector cable for EXC-02: CC02

\$25

This connector cable can be connected to EXC-02 and an Excel's



2P metal connector to modular adaptor: CC05

\$25

This adaptor can be connected to a detection sensor with a 2P metal connector and to a modular jack.



RCA cable to 2P metal connector conversion adaptor: CC08

\$25

This adaptor can be used to connect a RCA cable and 2P metal connector



2P metal connector to 3.5 mini plug adaptor: CC10

\$15

This adaptor can be connected to a detection sensor with a 2P metal connector and to a 3.5 mm mini plug.



3.5 mini plug to modular conversion adaptor: CC13

\$20

This adaptor can be connected to a detection sensor with 3.5 mini plug and to a modular jack.



2P metal connector to RCA cable adaptor: CC15

\$25

This adaptor can be used to connect a 2P metal connector and a RCA cable.



Modular jack to RCA cable conversion adaptor: CC17

\$15


This adaptor can be connected to a detection sensor with a modular jack and to a RCA cable



Hybrid Sensor cable for 2P metal connector: CC57-2

\$20

This cable can be connected to a hybrid detection sensor and to a 2P metal connector.



Adaptor Boxes for a Nurse Call System

The various caring sensors from Excel Engineering can be connected with other competitors' existing nurse call systems through the adaptor box.

*This product is only compatible to Excel Engineering sensors

* If you have an old nurse call system the adapter box may not be compatible. An additional fabrication fee may be applied.

Adapter Box for a wired connection



*NC-12-A

Adapter Box NC-10~


Type A \$55

Type B \$70

The adapter box can be turned off/on from the switch located on the top of the adapter.



*An optional bed hanger (BH-1) is available (the adaptor can be hung from the bed rail) BH-1



Adapter Box for a wireless connection Type A




*NC-5000NRX-12

Wireless Adapter Box NC-5000NRX-10~

Type A \$150

Type B \$160

This wireless adapter box has capabilities to receive detection signals and connect to an existing nurse call system. The wireless environment will help in creating a safe (no tripping over the cords) and more accessible caring environment.



Please refer to the following major nurse call system manufacturing companies' names and the type of nurse call connector pin to identify your specific nurse call connector pin. There are nurse call connector pins that may not be able to connect to Excel's adaptor box.

I. Nurse Call Connector Pins made by Aiphone

| Type: | A | A | A | A | A | A | A | B | B | B | B | A |
|----------------|------------|------------|---------|-------------------|--|--|--|------------------|-----------------|----------------------------|-------------------|--|
| Adaptor Box | NC-10 | NC-11 | NC-16 | NC-15 | NC-13-NEX | NC-13-NCX | NC-13-NBX | NC-14 | NC-17 | NC-18 | NC-19 | NC-22 |
| Product #: | 2P Regular | 3P Regular | IP Plug | AIPhone / Flat 2P | AIPhone / 4P NE model | AIPhone / 4P NC model | AIPhone / 4P NB model | AIPhone / Pla 5P | AIPhone / Din5P | AIPhone / for NFX option 3 | AIPhone / Din6P | AIPhone / Vi-nurse |
| Adaptor shapes | | | | | | | | | | | | |
| | | | | UAR-8 UBR-8 | * NE-RC, NE-RC/A: 5-E NCR-C: 5-C NB-RC/A: 5-B (there are some types that may not be compatible) | * NE-RC, NE-RC/A: 5-E NCR-C: 5-C NB-RC/A: 5-B (there are some types that may not be compatible) | * NE-RC, NE-RC/A: 5-E NCR-C: 5-C NB-RC/A: 5-B (there are some types that may not be compatible) | UA-SB8 UB-SB8 | NFR-8 | Option Plug | NF-SB NF-SB-TT | An optional input cord may be required (sold separately) |

2. Nurse Call Connector Pins made by CareCom

| Type: | A | A | A | A | A | B | A |
|----------------|------------|------------|---------|--|---|---------------------------------|-------------|
| Adaptor Box | NC-10 | NC-11 | NC-16 | NC-12-A | NC-12-D | NC-23 | NC-24 |
| Product #: | 2P Regular | 3P Regular | IP Plug | CareCom 4P analog | CareCom 4P digital | Well File: CareCom5P | CareComP/5P |
| Adaptor shapes | | | | | | | |
| | | | | No "Z" appears here. The color is ivory. BA-450F, BA-501 BA-440F, BA-511 | No "Z" appears here. The color is grey. BA-501Z BA-511Z | BA-601WZ BA-611WZ RB-830Z | |

*Please contact us if you cannot identify which 4P connection pin is yours.

*CareCom's "Well" series nurse call system requires a specific plug-in. The model number for the plug-in is BB-301WZ-2P/W6P-NC.

3. Nurse Call Connector Pins made by Panasonic


| Type: | A | A | B |
|----------------|---------------------------|---------|-----------------------|
| Adaptor Box | NC-20 | NC-16 | NC-21 |
| Product #: | Panasonic/ 2P grip button | IP Plug | Panasonic/ 2P handset |
| Adaptor shapes | | | |

If your nurse call connector pin is not listed, please contact us for consultation. There are nurse call connector pins that may not be able to connect to the Excel's adaptor box.

Adaptor Boxes for a Nurse Call System

SRH9000 product features: This is a wireless wandering/departure detection sensor system for individuals with dementia. When a care recipient's departure from the bed is detected the receiver starts an alarm or a melody sound to notify a caregiver. Since the system uses a 920 MHz radio transmission frequency with LoRa (Long Range Radio) technology, the wireless signal range could reach up to several kilometers in a straight line, and 100-300m within a building without a repeater.

SRH9000 with Base Receiver Package



When a care recipient steps on the mat an alarm/melody sound starts from the receiver to notify a caregiver.



When a care recipient gets up on the bed an alarm/melody sound starts from the receiver to notify a caregiver.

Special features: Use the sensor mat that serves as a bed and a floor sensor. A switch located on the transmitter can change the functions of the sensor either "getting up from the bed" or "leaving the bed" detection.




SRH9000 \$335


TAIS Code: 00682-000073

This package includes: One sensor mat, one transmitter, and one Base Receiver

SRH9000 with Portable Receiver Package



When a care recipient steps on the mat an alarm/melody sound starts from the receiver to notify a caregiver.



When a care recipient gets up on the bed an alarm/melody sound starts from the receiver to notify a caregiver.

Special features: Use the sensor mat that serves as a bed and a floor sensor. A switch located on the transmitter can change the functions of the sensor either "getting up from the bed" or "leaving the bed" detection.




SRH9000-P 325

TAIS Code: 00682-000072

This package includes: One Sensor mat, one transmitter, and one portable receiver

SH9000 with FloorSensor & Base Receiver Package



When a care recipient steps on the sensor an alarm/melody sound starts from the receiver to notify a caregiver.

SH9000 \$320

TAIS Code : 00682-000066

This package includes: One Floor sensor, one transmitter, and one base receiver



SR9000 with Bed Sensor & Base Receiver Package



When a care recipient gets up on the bed an alarm/melody sound starts from the receiver to notify a caregiver.


SR9000 \$305

TAIS Code : 00682-000065

This package includes: One Bed sensor, one transmitter, and one base receiver



SH9000 with Floor Sensor & Portable Receiver Package



When a care recipient steps on the sensor an alarm/melody sound starts from the receiver to notify a caregiver.


SH9000-P \$310

TAIS Code : 00682-000068

This package includes: One Floor sensor, one transmitter, and one portable receiver



SH9000 with Bed Sensor & Portable Receiver Package



When a care recipient gets up on the bed an alarm/melody sound starts from the receiver to notify a caregiver.


SR9000-P \$300

TAIS Code : 00682-000067

This package includes: One Bed sensor, one transmitter, and one portable receiver



WSR9000 with a Plate Sensor/Base Receiver Package



When a care recipient sits on the sensor an alarm/melody sound starts from the receiver to notify a caregiver.


SR9000-PS \$345

TAIS Code : 00682-000074

This package includes: One Plate sensor, one transmitter, and one base receiver



SR9000 with a Plate Sensor/Portable Receiver Package



When a care recipient sits on the sensor an alarm/melody sound starts from the receiver to notify a caregiver.

SR9000-PS-P \$335

TAIS Code : 00682-000071

This package includes: One Plate sensor, one transmitter, and one portable receiver



SH6900: This is a long-term care insurance rental approved product



See the instructional manual

SH6900 series is a long-term care insurance rental approved product.
SH6900 detects an individual's departure from the bed and departure from the room, and helps to prevent falling accidents.


SH6900 with Infrared Sensor to detect both "getting up from the bed" or "leaving the bed" activities.




- This wireless system's network coverage reaches up to approximately 100m in a straight line.
- The batteries are rechargeable. When the batteries need to be recharged the device informs a user.
- The Wandering/Departure Detection sensors can be registered as desired.
- The system is wireless, and no cables are required to install.

Specifications:

Wireless Transmission Method: 2.4 GHz Digital transmission Receiver (155g)
Dimensions: 125 x 90 x 23mm
Power: Lithium-Ion Polymer rechargeable battery built-in and an AC adapter included Transmitter (205g)
Dimensions: 98 x 89 x 85mm
Power: Lithium-Ion Polymer rechargeable battery built-in and an AC adapter included




Example: When a care recipient blocks the infrared...




An alarm sound starts from the receiver to notify a caregiver! The monitor displays the picture of the care recipient.

When a registered detection sensor in the transmitter is activated and the alarm sound starts on the receiver the receiver only displays the picture of the care recipient for 30 seconds. If the transmitter is not activated and no alarm sound starts the receiver does not display the picture of the care recipient.

SH6900 with Floor Sensor Package



When a care recipient steps on the mat...



An alarm sound starts from the receiver to notify a caregiver! The monitor displays the picture of the care recipient.

When a care recipient steps on the mat notify a caregiver. It is recommended to install the Floor sensor by the bed or the entrance door of the care recipient's room.

JAN Code : 4582107671100 TAIS Code : 00682-000035

SH6900 \$275




450 x 800 x 4.5 mm

SH6900 with Bed Sensor Package



When a care recipient gets up...




An alarm sound starts from the receiver to notify a caregiver! The monitor displays the picture of the care recipient.

When a care recipient gets up on the mat notify a caregiver. It is recommended to install the Bed sensor between a mattress pad and mattress just under the care recipient's shoulder blades when he/she lays down.


JAN Code : 4582107671117 TAIS Code : 00682-000036

SR6900 \$260




390 x 850 x 1.5 mm

SH6900 with Infrared Sensor Package



When a care recipient blocks the infrared...



An alarm sound starts from the receiver to notify a caregiver! The monitor displays the picture of the care recipient.




38 x 109 x 69 mm

When a care recipient blocks the infrared notify a caregiver. It is recommended to install the Infrared sensor by the bed or the entrance door of the care recipient's room.

JAN Code : 4582107671711 TAIS Code : 00682-000034


SH6900-MS2 \$295

SH6900 with Bed/Floor Sensor Package



When a care recipient gets up...

When a care recipient steps on the mat...



An alarm sound starts from the receiver to notify a caregiver! The monitor displays the picture of the care recipient.

This sensor mat serves as a bed and a floor sensor (i.e., "getting up from the bed" or "leaving the bed" detection).

JAN Code : 4582107671643 TAIS Code : 00682-000037

SRH6900-LS \$295



390 x 850 x 3.5 mm

Homecare Wireless Nurse Call System: NRM-6900



See the i nstructional




Transmitter/IP Camera

Receiver/Monitor

NRM-6900 \$170

JAN Code : 4582107671575



Call with your own voice

See the care recipient and can talk to her

Rechargeable

Voice and Live Video

Option available:
Long-lasting detection sensor

Option available:
Infrared Sensor

Option available:
Wandering prevention sensor


Option available:
Departure detection Sensor

- This is a wireless video monitoring system.
- The wireless signal range reaches approximately 100m in a straight line.
- *When a receiver leaves outside of the wireless signal reaching range the receiver displays a message and starts an alarm sound.
- Various detection sensors can be connected to this system wirelessly.
- The directions and angles of the transmitter/IP camera can be adjusted from the receiver/monitor.
- There is a thermal sensor built in.
- It is a 2-way wireless voice communication device.
- Transmitter/IP camera is rechargeable.

Time taken to be fully charged: approximately 3 hours for both receiver and transmitter when the devices are turned off.
Battery life when it is fully charged is approximately 4 hours for the receiver/monitor and 1.5 hours for the transmitter

Specifications:
Wireless Transmission Method: 2.4 GHz Digital transmission
Receiver (155g including a battery)
Dimensions: 125 x 90 x 23mm
Power: Lithium-Ion Polymer rechargeable battery built-in and an AC adapter included
Transmitter (205g including a battery)
Dimensions: 98 x 89 x 85mm
Power: Lithium-Ion Polymer rechargeable battery built-in and an AC adapter included

Detection Sensor Options



Homecare Wireless Nurse Call System HT-150



See the i nstructional

HT-150MRX \$50


*Transmitter is not included

JAN Code : 4582107671568

Additional transmitter
HT-150TX \$15

JAN Code : 4582107671384

Portable Call Button



Specifications
Wireless Transmission Method: Specified low power radio frequency Easy replacement battery - The transmitter uses a coin battery (CR2032).
Transmitter (weight: 25g)
Dimensions: 35 x 15 x 62mm
Power: Coin battery (CR2032) included
Receiver (weight: 129g)
Dimensions: 70 x 125 x 60mm
Power: A built-in battery and an AC adapter included

EXC-02S-B: On-site alert system for monitoring an individual departing from a bed or a wheelchair



- This is a long-term care insurance rental approved product
- This is an on-site alert device. It alerts a care recipient when he or she leaves from the detection sensor mat
- The detection sensor can be used in a wheelchair

EXC-02S-B \$200

TAIS Code : 00682-000083

Product Specifications

| Application | Call Request | | Monitor displays a care recipient and sensor detection | |
|---|---|--------------------------------|--|--|
| Name of Product | よへーる150 | | よへーる6900 | |
| | Multi Receivers | Transmitter | Receivers | Transmitter |
| Product Number | HT-150 | | NRM-6900 | |
| | HT-150MRX | HT-150TX | NRM-6900RX | NRM-6900TX |
| Price | | | | |
| Alert Method | — | Push button | — | Floor sensor Bed Sensor Infrared sensor |
| Alert Sound | Alarm Sound | — | Alarm Sound | — |
| Alert Sound Levels | 4 Levels | — | Fixed | — |
| Continuous Ringing Time | Approx. 20sec | — | Approx. 30sec | — |
| Signal Transmission Range (straight line) | Approx. 50m | | Approx. 100m | |
| Input | — | — | — | — |
| Power | AC adapter (DC5V/50mA) /designated rechargeable battery | One Coin battery CR2032 (3V) | AC adapter (DC5V/1A) /lithium ion polymer rechargeable battery | AC adapter (DC5V/1A) /lithium ion polymer rechargeable battery |
| Recharging Time | Approx. 8 hours | — | Approx. 3 hours | Approx. 3 hours |
| Battery Life | Approx. 30 hours | Approx. 2 years | Approx. 4 hours | Approx. 1.5 hours |
| Low Battery Indicator | Flashing the receiving light (every 5 sec) | Transmission light becomes dim | Displays a message on the screen | Power LED light flashes green |
| Accessories | — | — | One cable for monitor display (output) | |
| Dimensions (mm) | 70×60×125 | 35×15×62 | 125×90×23 | 98×89×85 |
| Weight (g) | 129 | 25 | 155 | 205 |
| Technical Number | — | 208-160024 | 208-170018 | |
| Materials | ABS | | | |

| Product Name | Cordless Infrared Sensor | Transmitter for sensors |
|--------------------------------|---|------------------------------|
| Product Number | RH-5000TXIR | RH-5000TXBF |
| Price | | |
| Alert Function | — | — |
| Alert Volume | — | — |
| Continuous Ringing Time | — | — |
| Display | — | Alarm light |
| Input | — | 3.5φ stereo jack |
| Receiving Channel | — | — |
| Wireless Specifications | Modulation method: GFSK Signal frequency bandwidth: 2.4GHz | |
| Power | AC adapter (DC9V/300mA) /9V alkaline battery | One Coin battery CR2032 (3V) |
| Battery Life | — | Approx. 1 year |
| Transmitter Battery Life Alert | Alarm notification | — |
| Accessories | — | — |
| Dimensions (mm) | 38×109×69 | 41×69×9 |
| Weight (g) | 140 | 23 |
| Materials | ABS | |

| Functions | Monitor displays a care recipient and sensor detection | | Sensor detection | | |
|---|--|--|-------------------------|---|---------------------|
| Product Name | わかーる6900 | | わかーる9000 | | |
| | Receiver | Transmitter | Base Receiver | Portable Receiver | Transmitter |
| Product Number | RH-6900 | | RH-9000RX | RH-9000PRX | RH-9000TXBF |
| | RH-6900RX | RH-6900TX | | | |
| Price | | | | | |
| Alert Method | — | Floor sensor Bed se nsor Infrared sensor | Sound Light | Sound Vibration | — |
| Alert Volume | Alarm sound | — | — | — | — |
| Alert Volume Adjustment | Fixed | — | 12 levels | Loud and quiet | — |
| Continuous Ringing Time | Approx. 30sec | — | Until reset | | — |
| Signal Transmission Range (straight line) | 約100m | | — | | |
| Input | — | — | — | — | 3.5φ stereo jack |
| Power | AC adapter (DC5V/1A) /lithium ion polymer rechargeable battery | AC adapter (DC5V/1A) /lithium ion polymer rechargeable battery | AC adapter (DC9V/300mA) | AC adapter (DC6V/300mA) /lithium ion polymer rechargeable battery | 2 AAA batteries |
| Recharging Time | Approx. 3 hours | Approx. 3 hours | — | — | — |
| Battery Life | Approx. 4 hours | Approx. 1.5 hours | — | — | Approx. 1 year |
| Low Battery Indicator | Display on the screen | Power LED light flashes green | Alarm light flashes | Display on the screen | Alarm light flashes |
| Dimensions (mm) | 125×90×23 | 98×89×85 | 160×140×30 | 60×91×20 | 50×85×18 |
| Weight (g) | 155 | 205 | 252 | 68 | 56 |
| Technical Number | 208-170015 | | — | | |
| 材 質 | ABS | | | | |

| System Name | NC5000 Wireless Nurse Call System | | | | | | | | |
|---------------------------------|---|---|--|----------------------|------------------------------|------------------------------|----------------------------|----------------------------|------------------------------|
| Name of Product | Base Receiver | Portable Receiver | Adapter Box | Call Button | Mini Transmitter | Door Sensor | Magnet Sensor Transmitter | Repeater | Mini Repeater |
| Product Number | NC-5000RX(C) | NC-5000PRX2 | NC-5000NRX-10~22 | NC-5000TXPB | NC-5000TXBF | NC-5000TXDR | NC-5000TXMG | NC-5000RP(C) | NC-5000RP2 |
| Price | | | | | | | | | |
| Alert Function | Sound & Light | Sound & Vibration | Nurse Call Connection | — | — | — | — | — | — |
| Alert Sound Volume | 12 levels | Loud & quiet | — | — | — | — | — | — | — |
| Continuous Ringing Time | Until reset | | — | — | — | — | — | — | — |
| Display | Power light Alarm light | Power light Monitor display | Power light 3 transmitter lights Pause light | Alarm light | | Alarm light | Transmitter light | Power light Alarm light | Alarm light |
| Input | — | — | — | Push button | 3.5φ stereo jack | — | — | — | — |
| Receiving Channels | 16 | 63 | 3 | — | — | — | — | — | — |
| Wireless Network Specifications | Modulation method: GFSK Signal frequency bandwidth: 2.4GHz | | | | | | | | |
| Power | AC adapter (DC9V/200mA) | AC adapter (DC6V/300mA) /lithium ion polymer rechargeable battery | AC adapter (DC5V/300mA) | 2 AAA batteries | One Coin battery CR2032 (3V) | One Coin battery CR2032 (3V) | 2 AAA batteries | AC adapter (DC9V/200mA) | Electrical outlet connection |
| Battery Life | — | — | — | Approx. 1 year | | Approx. 1 year | Approx. 1 year | — | — |
| Transmitter Battery Life Alert | Flashing alarm light | Display on the screen | Flashing transmitter light | Flashing alarm light | | Flashing alarm light | Flashing transmitter light | — | — |
| Accessories | Whip antenna | Belt clip strap | Bed hanger | — | — | Magnet | Hanger | Whip antenna | — |
| Dimensions | 160×140×30 | 60×91×20 | 82×157×33 | 50×85×18 | 41×69×9 | 41×69×9 | 62×31×111 | 160×140×30 | 77×86×37 |
| Weight | 252 | 68 | 200 | 56 | 23 | 23 | 140 | 240 | 90 |
| Materials | Housing:ABS | | | | | | | | |

| System Name | NC9000 Wireless Nurse Call System | | | | |
|---------------------------------|-------------------------------------|---|----------------------------|----------------------|----------------------------|
| Product Name | Base Receiver | Portable Receiver | Adapter Box | Transmitter | Magnet Sensor Transmitter |
| Product Number | RH-9000RX | RH-9000PRX | NC-9000GW | RH-9000TXBF | NC-9000TXMG |
| Price | | | | | |
| Alert Function | Sound & light | Sound & vibration | — | — | — |
| Alert Sound Volume | 10 levels | Loud & quiet | — | — | — |
| Continuous Ringing Time | Until reset | | — | — | — |
| Display | Power light Alarm light | Charging light & display | Power light Alarm light | Alarm light | Transmitter light |
| Input | — | — | — | 3.5φ stereo jack | — |
| Receiving Channels | 16 | 100 | — | — | — |
| Wireless Network Specifications | Signal frequency bandwidth: 9200MHz | | | | |
| Power | AC adapter (DC5V/300mA) | AC adapter (DC6V/300mA) /lithium ion polymer rechargeable battery | AC adapter (DC5V/300mA) | 2 AAA batteries | 2 AA batteries |
| Battery Life | — | — | — | Approx. 1 year | |
| Transmitter Battery Life Alert | Flashing alarm light | Display on the screen | — | Flashing alarm light | Flashing transmitter light |
| Accessories | — | Belt clip strap | — | — | Hanger |
| Dimensions (mm) | 160×140×30 | 60×91×20 | 160×140×30 | 50×85×18 | 62×31×111 |
| Weight (g) | 252 | 68 | 240 | 56 | 140 |
| Materials | Housing:ABS | | | | |

| System Name | NC5000 | NC9000 |
|---------------------------------|---|------------------------------|
| Product Name | Nurse call button | Diaper sensor transmitter |
| Product Number | NC-5000TXNCB | NC-5000TXWA |
| Price | | |
| Input | Push button | Metal clump |
| Wireless Network Specifications | Modulation method: GFSK Signal frequency bandwidth: 2.4GHz | |
| Battery Life | Approx. 1 year | Approx. 1 year |
| Power | One Coin battery CR2032 (3V) | One Coin battery CR2032 (3V) |
| Accessories | Strap | — |
| Dimensions (mm) | 112×50×30 | 43×53×13 |
| Weight (g) | 44 | 29 |
| Materials | Housing: ABS | Housing: ABS |

| Product Name | Magnet sensor |
|-------------------------|---|
| Product Number | EMG-01 |
| Price | |
| Alert Function | Sound |
| Alert Volume | High: approx. 95dBA Low: approx. 85dBA |
| Continuous Ringing Time | Until magnet reattaches |
| Display | Battery life display light |
| Power | 3 AA batteries |
| Battery Life | Approx. 1 year |
| Low Battery Display | Flashing battery life display light |
| Accessories | Hanger |
| Dimensions (mm) | 62×31×111 |
| Weight (g) | 160 |
| Materials | Housing: ABS |

About Our Company

About Our Company

Name of Company: Excel Engineering Co., LTD

Address: Seika Building 9F, 2-12-14 Ogawamachi, Chiyoda-ku, Tokyo, Japan 101-0052

Phone: 03-5280-7120 **Fax:** 03-5280-7123

Date Founded: December 1996

Capital : 50 million yen

President and CEO: Katsunori Yasuoka

The banks that we do business with: Mitsubishi-TokyoUFI Bank Shintomicho branch, Mizuho Bank Kanda

Ogawacho Branch, and Mitsui-Sumitomo Bank Kanda Branch

History

December, 1996

Excel Engineering Co., LTD was founded

February, 1997

Excel Engineering Hong Kong was founded

April, 1999

"Leaving from the bed" Detection Sensor was first sold

April, 2000

Act for Enforcement of the Long-Term Care Insurance Act

April, 2002

The Home Care Wireless Monitoring Systems NRM series and Wireless Diaper Sensor (WS-900) were first sold

May, 2002

"Departing from the bed" Twin-pad detection sensor package was first sold

June, 2003

Excel Engineering Osaka Branch was opened NRM190 was first sold

March, 2003

Excel Engineering main office was moved to Shintomi, Chuo-ku, Tokyo

October, 2004

Pad-Less Wheelchair Detection Sensor was first sold Wireless Nurse Call System was first sold

May, 2005

Excel Engineering main office was moved to Muromachi, Nihonbashi, Chuo-ku, Tokyo

October, 2005

Takamatsu branch was opened

October, 2006

Infrared Double Sensor was first sold

December, 2008

Wireless Nurse Call System (NC-5000) was first sold

April, 2009

Nagano Distribution Center was opened

January, 2014

Takamatsu office was relocated due to business expansion

October, 2014

Exhibited at the International Home Care & Rehabilitation Exhibition (H.C.R) for the 10th consecutive year

November, 2014

Exhibited at MEDICA trade Fair in Germany for the 5th consecutive year

March, 2015

Started repair service for care facility products at Nagano office

April, 2015

Established Yobe-ru CO., LTD which provided home care products

May, 2015

Contracted with Bio Silver CO., LTD to expand sales business for products related to "getting up from the bed" and "departing from the bed" detection system

April, 2016

Yobe-ru CO., LTD' s home care product was first sold

October, 2016

Yobe-ru CO., LTD, Tokyo office was open

February, 2017

Contract with Bio Silver CO., LTD ended

May, 2017

Excel Engineering main office was moved to Kanda-Ogawamachi, Chiyoda-ku, Tokyo

April, 2018

Cloud-Based Total Care System was first sold

October,2018

NC9000 infrared detection sensor system was first sold

Company' s Information

Excel Engineering Co., LTD.

<http://www.excel-jpn.com>

Headquarter/Tokyo Office: Seika Building 9F, 2-12-14 Ogawamachi, Chiyoda-ku, Tokyo, Japan 101-0052

Phone: 03-5280-7120 **Fax:** 03-5280-7123

Takamatsu Office: 681-1 Haze-Chou, Takamatsu, Kagawa, Japan 761-8054

Phone: 087-815-0805 **Fax:** 087-869-2701

Nagano Office (Facility Repair Department): 495-36 Akiwa-Toiya-Machi, Ueda, Nagano, Japan 386-0041

Phone: 0268-75-7140 **Fax:** 0268-75-7141

Yobe-ru Co., LTD. (Product Sale Department)

<http://www.yobe-ru.jp>

Headquarter/Tokyo Office: Seika Building 9F, 2-12-14 Ogawamachi, Chiyoda-ku, Tokyo, Japan 101-0052

Phone: 03-5280-7120 **Fax:** 03-5280-7123

Sales Division: 681-1 Haze-Chou, Takamatsu, Kagawa, Japan 761-8054

Phone: 087-815-0805 **Fax:** 087-869-2701

Distribution Center: 495-36 Akiwa-Toiya-Machi, Ueda, Nagano, Japan 386-0041

Phone: 0268-75-7140 **Fax:** 0268-75-7141

| Product Name | Hybrid Floor Sensor Size: Small | Hybrid Floor Sensor Size: Medium | Hybrid Floor Sensor Size: Large | Cordless Bed Sensor Size: Wide | Bed Sensor Size: Wide | Bed Sensor Size: Regular |
|-------------------|--|-------------------------------------|------------------------------------|-----------------------------------|--------------------------|-----------------------------|
| Product Number | EFM-57/S-C | EFM-57/M-C | EFM-57/L-C | EPB58-W-2 | EPB54-W-I | EPB54 |
| Price | | | | | | |
| Output | 3.5φ mono plug / 2P metal connector / 4/4 modular plug | | | 3.5φ mono plug | 4/4 modular plug | |
| Dimensions | 450×800×4.5 | 600×900×4.5 | 600×1200×4.5 | 390×850×1.5 | 390×767×6.5 | 125×765×1.5 |
| Weight | 880 | 1550 | 2200 | 500 | 540 | 200 |
| Materials | PVC | | | | | |
| Cable Length (mm) | 2500 | | | — | 2000 | |

| Product Name | Small Cordless Side-Rail Sensor Pad | Cordless Side- Rail Sensor Pad | Cordless Long-Life Detection Mat | Long-Life Detection Mat Wheel | Wheel Chair Sensor Pad | Toilet Sensor |
|-------------------|--|-----------------------------------|--|-------------------------------------|---------------------------|----------------|
| Product Number | EPB58-SR/S | EPB58-SR/L | EBF-55 | EBF-51/S-I | EPC53 | ETP56-5 |
| Price | | | | | | |
| Output | 3.5φ mono plug | | | 4/4 modular plug | | 3.5φ mono plug |
| Dimensions | 170×280×8 | 260×700×8 | 390×850×3.5 | 381×765×3.5 | 354×382×6.5 | 82×38×2 |
| Weight | 148 | 490 | 1200 | | 220 | 28 |
| Materials | PVC | | | | | |
| Cable Length (mm) | —— | | | 2000 | 750 | 800 |

| | | | |
|-------------------------|---|---|---|
| Product Name | Infrared Sensor | Infrared Sensor MS3000 Transmitter built in | Infrared Sensor MS9000 Transmitter built in |
| Product Number | MS-3100D | MS-3100E | NC-9000TXMS |
| Price | | | |
| Alert Function | Nurse call | | |
| Display | LED lights: Power, Detection, Alert | | |
| Detection Capacity | Area sensor: approx. 5m, 30 degrees Infrared sensors: approx. 1.3~1.75m, 3 beams | | |
| Wireless Specifications | _____ | Modulation method: GFSK Signal frequency bandwidth: 2.4GHz | Signal frequency bandwidth: 9200MHz |
| Output | Make contact | | |
| Power | AC adapter (DC5V/300mA) | | |
| Accessories | NC connection cord | _____ | |
| Dimensions | 155×110×40 | | |
| Weight | 300 | 310 | 310 |
| Materials | Housing: ABS | | |

| Product Name | Wired Nurse Call Activation/ In-room Alert System |
|-------------------------|---|
| Product Number | EXC-02B |
| Price | |
| Alert Function | Sound, light, nurse call activation |
| Alert Volume | Loud: approx. 100dBA Medium: approx. 85dBA |
| Continuous Ringing Time | Until reset |
| Display | Power light, Alarm |
| Input | 4/4 modular plug |
| Low Battery Indicator | 6φ mono jack Contact pressure resist: 60V Contact current: 400mA ON resistance: 7Ω |
| Power | AC adapter (DC6V/100mA) /2 AAA batteries |
| Battery Life | Approx. 1 year |
| Low Battery Indicator | Flashing alarm light |
| Accessories | Bed hanger |
| Dimensions (mm) | 60×25×120 |
| Weight (g) | 113 |
| Materials | Housing: ABS |

| Functions | Sensor Accessories | | | | | |
|-------------------|-----------------------------------|---------------------|---------------------|----------------|----------------------------------|-----------------------|
| Product Name | Nurse call connection cable | Long stand | Short stand | Bed holder | Adjustable angle connector | 2PNC cable |
| Product Number | CC02 | MS-3000-2 | MS-3000-3 | BH-2 | PT-I | MS-03 |
| Price | | | | | | |
| Output | 2P metal connector | — | — | — | — | 2P metal connector |
| Dimensions (mm) | — | Height: 600~1100 | Height: 600~1100 | 100×170 ×70 | 41×68 | — |
| Weight (g) | 70 | 2700 | 1400 | 165 | 68 | 90 |
| Materials | — | Fe | Fe | ABS | — | — |
| Cable Length (mm) | 2000 | — | — | — | — | 2000 |



MEMO

Handwriting practice area on page 37, consisting of 20 horizontal dashed lines.



MEMO

Handwriting practice area on page 38, consisting of 20 horizontal dashed lines.